


# BIOHORIZONS®

## Pro surgical kit (PRO5000)

Includes the instruments required to place BioHorizons Tapered Implants.†

### 1. Guide Fixation



Pilot Drill	Fixation Pin
122-024	CGS-FP

### 2. Pre Drills



Ø 2.8mm	Ø 3.2mm	Ø 3.7mm	Ø 4.1mm	Ø 4.7mm	Ø 5.4mm
PROPRE25	PROPRE28	PROPRE32	PROPRE37	PROPRE41	PROPRE47

### 3. Ancillary Instruments



050" (1.25mm) Hex Driver	4mm Square Drive Extender	Drill Extender	1.5mm starter drill	Tapered HD Drill, 2.0mm
135-351	300-205	122-100	122-015	TSD2020HD

### 4. Parallel Pins



Straight Parallel Pins	20° Angled Parallel Pins
144-100	144-200

### 5. Bone Levelers



Grey	Yellow	Purple
PROBLGY	PROBLYL	PROBLPR

### 6. Guide Adapters



Grey	Yellow	Purple
PROGAGY	PROGAYL	PROGAPR

### 7. Square Converter



4mm Square Drive Converter
CGS-4SC

### 8. Implant Drivers




Grey	Yellow	Purple	Grey	Yellow	Purple
PROIDGY	PROIDYL	PROIDPR	PROCDGY	PROCDYL	PROCDPR

### 9. Tissue Punches




Grey	Yellow	Purple
PROTPGY	PROTPYL	PROTPPR



Freehand = width-increasing drills

	Soft Bone	Dense Bone
Ø3.0	—	Ø2.5
Ø3.3/3.4	—	Ø2.8
Ø3.8	Ø2.8	Ø3.2
Ø4.2	Ø3.2	Ø3.7
Ø4.6	Ø3.7	Ø4.1
Ø5.2	Ø4.1	Ø4.7
Ø5.8	Ø4.7	Ø5.4

### 13. Conical Implant Drivers



Narrow, Handpiece	Narrow, Ratchet
CNIDH	CNIDR
Regular, Handpiece	Regular, Ratchet
CRIDH	CRIDR

### 10. Length Drills



	2.5mm	2.8mm	3.2mm	3.7mm	4.1mm	4.7mm	5.4mm
9mm	PROLD2509	PROLD2809	PROLD3209	PROLD3709	PROLD4109	PROLD4709	PROLD5409
10.5mm	PROLD2510	PROLD2810	PROLD3210	PROLD3710	PROLD4110	PROLD4710	PROLD5410
12mm	PROLD2512	PROLD2812	PROLD3212	PROLD3712	PROLD4112	PROLD4712	PROLD5412
15mm	PROLD2515	PROLD2815	PROLD3215	PROLD3715	PROLD4115	PROLD4715	PROLD5415

### 11. Extended Shank Drills




2.5mm	2.8mm	3.2mm	3.7mm	4.1mm	4.7mm	5.4mm
TSD4025HD	TSD4028HD	TSD4032HD	TSD4037HD	TSD4041HD	TSD4047HD	TSD4054HD

### 14. Internal Hex Implant Drivers



3.0mm Handpiece, Regular	3.0mm Ratchet, Regular	3.5/4.5mm Handpiece	3.5/4.5mm Ratchet
TP3IDHR	TP3IDRR	TYGIDH	TYGIDR

### 12. Ratchet



Ratchet	130-000
---------	---------

customer care: 888.246.8338  
shop online at [store.biohorizons.com](http://store.biohorizons.com)

†Mount-free Tapered Pro Conical, Tapered Pro, Tapered Internal, Tapered Plus, Tapered 3.0. Spare slots can be used to hold preferred surgeon instrumentation. Refer to L02087 for further information and instructions. ©BioHorizons. All Rights Reserved. Images shown not to scale. L02088 Rev B MAY 2024